

Exp #8  
CHM 1045

Important info

**Part A**

1. Weigh 1 g unknown, add 150 mL DI H<sub>2</sub>O, stir 3 minutes (DO NOT adjust pH of water)
2. Heat (but DO NOT BOIL!)
3. We will use gravity filtration
4. Make sure sol'n is clear before extracting some
5. Ok
6. Dry in oven for 20 min

**Part B**

1. I'll show you how to use the centrifuge
2. Ok
3. Ok

Pre-lab Questions To Do (remember to have them done before lab!)

1, 4

Post-lab Questions To Do

2, 3,5